

# RAHM Sensor Box Instructions

## Required Materials

- RAHM box
- Cell Guardian device
- Laptop or other computer with wireless network support

## General Setup

- 1) Connect the Cell Guardian to the RAHM box using the Gray (Device) cable
- 2) Turn on the computer in the RAHM box
  - a. It will take up to 5 minutes for the system to come online
- 3) Connect to the RAHM box wireless network with the laptop
  - a. SSID: RAHM\_DEMO
  - b. Password: Guardian2025
- 4) Open a browser on the laptop
- 5) Go to <http://192.168.42.1>
  - a. This should present a directory of web pages
- 6) Click on the Cell Guardian interface link
- 7) At the Keycloak login page enter the following credentials
  - a. Username: admin
  - b. Password: Guardian2025
- 8) After logging in you should arrive at the Cell Guardian interface
- 9) Click on the active device

## Shutting Down

The computer should be turned off before removing power.

- 1) Short press the power button
- 2) Wait for the power LED to turn off
- 3) Unplug the RAHM box

## Connecting to the Internet

The red cable on the RAHM box provides internet connectivity. Plug the cable into a router, hotspot, or other DHCP-enabled ethernet port. The RAHM box is remotely accessible when the red cable is connected to the Internet.